



## Bonus Activity #1: Pipette by Number

Collect a piece of Parafilm, wax paper, plastic sheet protector, or Petri dish. Using the P20, P200, and P1000 pipettes, recreate the designs below using water dyed with food coloring.



We would love to see your artistic creations! Send images of pipette art to [info@wolbachiaproject.com](mailto:info@wolbachiaproject.com) to be featured on our social media accounts, website, and/or newsletter.

Would you like to create a new template for this activity? Indicate the droplet sizes and we will add your design(s) to this lab guide.





## Bonus Activity #2: Freestyle Pipette Art

### PREPARATION

1. Print out one or more of the following template pages and cut into quarters. Each quarter (1/4) page should feature a single art template.
2. Cut a plastic sheet protector into quarters. This should yield eight small 1/4's. Alternatively, you may use Petri dishes.
3. Prepare small cups of water dyed with food coloring.

### ACTIVITY

1. Collect an art template and place it under a plastic sheet protector or Petri dish.
2. Using the P20, P200, and P1000 pipettes, recreate the design using colored water. Change droplet sizes to enhance the pipette art.

